



WELL life™

Strep A Test

Group A Streptococci (GAS) are one of the most common bacteria encountered daily, and they result in acute infections with a wide array of manifestations in both adults and children. Our rapid, lateral-flow diagnostic tool is useful for physicians to identify the presence of GAS in an office setting and to start the appropriate treatment.

- Results available at 10 minutes
- 94.7% sensitivity;
98.5% specificity vs. culture
- Room temperature storage
- CLIA Waived
- External controls included



Item Code	Description	Box Size / Case Size	Complexity
PRE-STa-25	25-Strip, Step A Throat Swab Test	25 / 20	CLIA WAIVED



WELL life™

iFOB

Immunochemical Fecal Occult Blood Test (FIT)

Fecal Immunoassay Tests (FIT) are also commonly referred to as Immunoassay Fecal Occult Blood Tests (iFOB) and are recommended to replace the older, guaiac-based stool cards (gFOB) for the detection of lower gastrointestinal (G.I.) bleeding. Physicians often use FIT, along with other indicators, as a method of Colorectal Cancer screening for average risk populations, and other lower G.I. disorders.

- Results available at 10 minutes
- Analytical sensitivity 45 ng/mL
- Specific to human hemoglobin
- Room temperature storage
- CLIA Waived



Item Code	Description	Box Size / Case Size	Complexity
PRE-FIT-25KIT	25 Test cassettes, 25 Collection tubes, 25 Mailers	25 / 18	CLIA WAIVED
PRE-FIT-25A	25 Test cassettes, 25 Collection tubes	25 / 20	CLIA WAIVED

* Ordered upon request

Rel: 2023/05/09