





WELL life™

# iFOB

## Immunochemical Fecal Occult Blood Test (FIT)

Fecal Immunoassay Tests (FIT) are also commonly referred to as Immunoassay Fecal Occult Blood Tests (iFOB) and are recommended to replace the older, guaiac-based stool cards (gFOB) for the detection of lower gastrointestinal (G.I.) bleeding. Physicians often use FIT, along with other indicators, as a method of Colorectal Cancer screening for average risk populations, and other lower G.I. disorders.

- Results available at 10 minutes
- Analytical sensitivity 45 ng/mL
- Specific to human hemoglobin
- Room temperature storage
- CLIA Waived



Item Code	Description	Box Size / Case Size	Complexity
PRE-FIT-25KIT	25 Test cassettes, 25 Collection tubes, 25 Mailers	25 / 18	CLIA WAIVED
PRE-FIT-25A	25 Test cassettes, 25 Collection tubes	25 / 20	CLIA WAIVED

\* Ordered upon request

Rel.: 2023/05/09