



INDICAID™

INDICAID™ COVID-19 Rapid Antigen Test



The INDICAID™ COVID-19 Rapid Antigen Test is a point-of-care kit designed for the detection of SARS-CoV-2 antigens in direct nasal swab samples.

The test is intended for use by healthcare professionals in symptomatic individuals within 5 days after onset of symptoms.


This product has received FDA Emergency Use Authorization.

Product Advantages

- **Fast:** Results in 20 minutes
- **Convenient:** No equipment or training needed
- **Simple Sampling:** Self-collected shallow nasal samples
- **Batch Testing:** Easily collect and test multiple samples


Intuitive Workflow

01




Collect Nasal Swab

02



Insert and Twist

03



Apply Sample


04



Read Visually

Batch Testing

Intuitive design allows for the concurrent collection of a large number of samples followed by batch testing of individual samples within 2 hours.



Collection

Testing

Clinical Performance

In a prospective US clinical study, INDICAID™ accurately identified 84% of those who were positive (PPA) and 97% of those who were negative (NPA) for SARS-CoV-2.

Product performance against variants of concern are evaluated on an ongoing basis.

Technical Bulletin available upon request.

About Us

PHASE Scientific is a high-growth biotech company founded by bioengineers from UCLA. We build tools that empower people by giving them better information about their health.

indicaid@phasesci.com

us.phasescientific.com

PHASESCIENTIFIC