



## HbA1c results in minutes

Informed conversations at the point of care

The A1CNow<sup>®</sup> system from PTS Diagnostics provides healthcare professionals with a fast and easy way of getting HbA1c results—allowing informed conversations about how patients are managing their diabetes in minutes, not days.

For more information, visit [www.a1cnow.com](http://www.a1cnow.com).

### Key Features


- National Glycohemoglobin Standardization Program certified
- CLIA-waived
- Tests may be reimbursable

# HbA1c results in minutes

Informed conversations at the point of care



A1CNow<sup>®</sup>+ system  
(10-count test kit)  
3024



A1CNow<sup>®</sup>+ system  
(20-count test kit)  
3021

## A1CNow<sup>®</sup>+

### More efficient than the lab

- Results in 5 minutes
- Small (5 $\mu$ L) blood sample from fingerstick
- Portable – use in multiple exam rooms
- Easy to use – minimal training required
- Tests may be reimbursable



Website: [ptsdiagnostics.com](http://ptsdiagnostics.com)  
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608  
Toll-Free in the US: 1 (877) 870-5610