

Lipid and Glucose Results at the Point of Care



On-site lipid and glucose results
in as little as 90 seconds

Compact and lightweight,
easy to transport and store

Handheld, battery-operated

Room temperature test strips

Fingerstick with small sample

Easy data transfer

Rapid throughput and low waste,
with minimal start-up cost

Fully reimbursable



Simultaneous testing
of lipid profile and glucose

CLIA-waived
Meets NCEP guidelines for accuracy and precision
FCC and ISO certified wireless communications
CRMLN-certified for CHOL and HDL tests
FDA-cleared

The CardioChek Plus Lipid+eGLU® Smart Bundle Pack™ Includes:

PTS Panels® lipid panel test strips and eGLU test strips



Lipid Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Triglycerides

Calculates: LDL Cholesterol, Total Cholesterol/HDL Ratio, Non-HDL Cholesterol



CHOL+HDL+GLU Panel Test Strip

Measures: Total Cholesterol, HDL Cholesterol & Glucose

Calculates: Total Cholesterol/HDL Cholesterol Ratio, Non-HDL Cholesterol



METABOLIC CHEMISTRY Panel Test Strip

Measures: HDL Cholesterol, Triglycerides & Glucose



CHOL+GLU Panel Test Strip

Measures: Total Cholesterol & Glucose



CHOL+HDL Panel Test Strip

Measures: Total Cholesterol & HDL Cholesterol

Calculates: Total Cholesterol/HDL Cholesterol Ratio



TOTAL CHOLESTEROL Test Strip

Measures: Total Cholesterol



HDL CHOLESTEROL Test Strip

Measures: HDL Cholesterol



TRIGLYCERIDES Test Strip

Measures: Triglycerides



GLUCOSE Test Strip

Measures: Glucose



eGLU Test Strip

Measures: Electrochemical Glucose

Single-analyte test strips available as special order include Creatinine and Ketone.



PTS Diagnostics
7736 Zionsville Road, Indianapolis, IN 46268 USA
Office: +1 (317) 870-5610 FAX: +1 (317) 870-5608
Toll-Free in the US: 1 (877) 870-5610